

## The new Series 7000 Large Power System from MicroAire®

Offering the ultimate in convenience and versatility, Series 7000 hand pieces function as both cordless and electric instruments. Choose between the industry's first autoclavable or aseptic transfer NiMH battery or the Direct Power Console.

### MicroAire BATTERY CHARGER 7500-700 MicroAire LARGE BATTERY PACK 7505-710

Features **14.4 Volt Nickel Metal Hydride (NiMH) rechargeable, autoclavable or aseptic transfer battery**

MicroAire NiMH batteries have 30% more power than the leading NiCad system

Four charging bays

NiMH batteries do not contain memory and have longer run times

Charges all MicroAire batteries within 90 minutes

MicroAire NiMH batteries last 3x as long as the competition

### DIRECT POWER CONSOLE 7025

Features **Simple dedicated electric power supply converts hand pieces to electric**

Fully autoclavable  
DPC Cable

Use as dedicated electric or as a battery backup

### ASEPTIC BATTERY SYSTEM

Features **Non-sterile aseptic transfer, battery never needs autoclaving**

For use on any MicroAire battery powered hand piece

## Instruments & Accessories



Catalog Number	Description
7025	Direct Power Console (DPC)
7007-000	DPC Console Pack
7100-001	Zimmer® Reamer Coupler (Hudson® Modified Trinkle), 250 RPM
7100-002	Trinkle Coupler, 750 RPM
7100-003	Jacobs Coupler 1/4", 750 RPM
7100-004	Keyless Coupler 1/4", 750 RPM
7100-005	Keyless Coupler 1/8", 750 RPM
7100-006	Synthes® Quick Connect Coupler, 750 RPM
7100-008	Hudson® Reamer Coupler, 250 RPM
7100-010	Jacobs Coupler 1/4", 250 RPM
7100-011	Harris/Aesculap® Reamer, 250 RPM
7100-111	Synthes® Reamer Coupler, 250 RPM
7100-112	Synthes® Non-Reversing Reamer Coupler, 250 RPM
7505-045	Wire/Pin Driver Coupler, 1.6mm- 2.8mm
7505-050	Pin Driver Coupler, 3.2mm- 4.0mm
7100-051	Zimmer® Coupler (Hudson® Modified Trinkle), 750 RPM
7100-058	Hudson® Coupler, 750 RPM
7100-474	Pulse Lavage Coupler
7505	Drill/Reamer
7500-700	Universal Battery Charger
7505-710	Large Power Battery Pack
7506	Oscillating Saw
7507	Reciprocating Saw
7505-040	Sterilization Case
7500-615	Aseptic Battery Housing
7500-620	Aseptic Battery Pack
7500-625	Charging Adapter
7500-630	Aseptic Transfer Shield

See MicroAire Knee Blade brochure for a complete list of cutting accessories.

## MICROAIRE®

MicroAire Surgical Instruments LLC  
1641 Edlich Drive  
Charlottesville, Virginia 22911  
U.S.A.

(800) 722 0822  
(434) 975 8000  
Fax: (800) 648 4309  
Fax: (434) 975 4131  
www.microaire.com



Product Information

SERIES 7000™

MICROAIRE®



**7505 Drill/Reamer**

**DRILL**

Speed **750 RPM**  
Torque **50 in.lbs (5.6 Nm)**

**REAMER**

Speed **250 RPM**  
Torque **150 in.lbs (16.9 Nm)**



**7506 Oscillating Saw**

**SAW**

Speed **15,000 CPM**  
Arc **4°**  
8 Position Rotating Head  
Two Speed Control



**7507 Reciprocating Saw**

**SAW**

Speed **15,000 CPM**  
Blade  
Stroke **3.2mm**  
360° Blade Rotation

**Drilling**



**Keyless Coupler 1/4"**  
Speed **750 RPM**  
Catalog # **7100-004**



**Keyless Coupler 1/8"**  
Speed **750 RPM**  
Catalog # **7100-005**



**Synthes® Quick Connect Coupler**  
Speed **750 RPM**  
Catalog # **7100-006**



**Zimmer® (Hudson®) Modified Trinkle Coupler**  
Speed **750 RPM**  
Catalog # **7100-051**



**Hudson® Coupler**  
Speed **750 RPM**  
Catalog # **7100-058**



**Trinkle Coupler**  
Speed **750 RPM**  
Catalog # **7100-002**

**Reaming**



**Zimmer® (Hudson®) Modified Trinkle Reamer Coupler**  
Speed **250 RPM**  
Catalog # **7100-001**



**Hudson® Reamer Coupler**  
Speed **250 RPM**  
Catalog # **7100-008**



**Jacobs® Coupler 1/4"**  
Speed **250 RPM**  
Catalog # **7100-010**



**Harris/Aesculap® Reamer Coupler**  
Speed **250 RPM**  
Catalog # **7100-011**



**Synthes® Reamer Coupler**  
Speed **250 RPM**  
Catalog # **7100-111**



**Synthes® Non-Reversing Reamer Coupler**  
Speed **250 RPM**  
Catalog # **7100-112**

**Pulse Lavage**



**Pulse Lavage Coupler**  
Catalog # **7100-474**  
Uses 5740 tubing sets

**Pin Driving**




**Jacobs® Coupler 1/4"**  
Speed **750 RPM**  
Catalog # **7100-003**



**Wire/Pin Driver Coupler**  
Speed **750 RPM**  
Catalog # **7505-045**



**Pin Driver Coupler**  
Speed **750 RPM**  
Catalog # **7505-050**



**Precision, Quality, MicroAire**